

## Product Datasheet

# Feline IL-12 p40 Recombinant Protein (orb1972929)

<b>Description</b>	The Feline IL-12/IL-23 p40 yeast-derived recombinant protein is not tagged and is naturally endotoxi...
<b>Target</b>	IL-12 p40
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	IWELEKNVYV VELDWHPDAP GEMVLT CNT PEEDDITWTS DQSSEVLGSG KTLTIQVKEF ADAGQYTCHK GGEVLSHSFL LIHKKEDGIW STDILREQKE SKNKIFLKCE AKNYSGRFTC WWLTAISTDL KFTVKSSRGS SDPQEVTCGA ATLSAEKVRV DNRDYKKYTV ECQEGSACPA AEESLPIEVV VDAIHKLKYE NYTSSFFIRD IIKPDPPKNL QLKPLKNSRH VEVSWEYPDT WSTPHSYFSL TFGVQVQGN NREKKDRLSV DKTSAKVVCH KDAKIRVQAR DRYSSSSWSN WASVSCS (307)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	35.0 kDa
<b>Entrez</b>	<b>768273</b>
<b>Expiration Date</b>	6 months from date of receipt.