

Product Datasheet

Bovine Prolactin Receptor Recombinant Protein (orb1972932)

Description	The Bovine Prolactin Receptor yeast-derived recombinant protein is not tagged and is naturally endot...
Target	PRLR
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QSPPEKPKLV KCRSPGKETF TCWWEPEGADG GLPTNYTLTY HKEGETLIHE CPDYKTGGPN SCYFSKKHTS IWKMYVITVN AINQMGISSS DPLYVHVITYI VEPEPPANLT LELKHPEDRK PYLWIKWSPP TMTDVKSGWF IIQYEIRLKP EKATDWETHF TLKQTQLKIF NLYPGQKYL QIRCKPDHGY WSEWSPESSI QIPNDFPVKD TSM (213)
Purity	98%
Source	Yeast
MW	24.8 kDa
Entrez	281422
Expiration Date	6 months from date of receipt.