

## Product Datasheet

# Bovine Prolactin Receptor Recombinant Protein (orb1972933)

<b>Description</b>	The Bovine Prolactin Receptor yeast-derived recombinant protein is not tagged and is naturally endot...
<b>Target</b>	PRLR
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	QSPPEKPKLV KCRSPGKETF TCWWEPEGADG GLPTNYTLTY HKEGETLIHE CPDYKTGGPN SCYFSKKHTS IWKMYVITVN AINQMGISSS DPLYVHVITYI VEPEPPANLT LELKHPEDRK PYLWIKWSP TMTDVKSGWF IIQYEIRLKP EKATDWETHF TLKQTQLKIF NLYPGQKYL QIRCKPDHGY WSEWSPESSI QIPNDFPVKD TSM (213)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	24.8 kDa
<b>Entrez</b>	<a href="#">281422</a>
<b>Expiration Date</b>	6 months from date of receipt.