

Product Datasheet

Chicken IL-34 Recombinant Protein (orb1972936)

Description	The Chicken IL-34 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-34
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ECELARLLQD KLR YEMRLQY MKHNFIDYT LRVQHEEVL R TANVTRLRDG KVSEASLRYL WFHACSQAVL HILEVLPEKH PSRGYTQELS QLLDALGVEY SGYRQSDVDA VVADLVKQLH SGDSRQKAVR PKALLDNCLK VLRMLFGAHC RWDSA (155)
Purity	98%
Source	Yeast
MW	18.0 kDa
Entrez	100858424
Expiration Date	6 months from date of receipt.