

## Product Datasheet

# Mouse Prolactin Recombinant Protein (orb1972937)

<b>Description</b>	The Mouse Prolactin yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
<b>Target</b>	PRL
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	LPICSAGDCQ TSLRELFDRV VILSHYIHTL YTDMFIEFDK QYVQDREFMV KVINDCPTSS LATPEDKEQA LKVPPEVLLN LILSLVQSSS DPLFQLITGV GGIQEAPPEYI LSRAKEIEEQ NKQLLEGVEK IISQAYPEAK NGGIYFVWSQ LPSLQGVDEE SKILSLRNTI RCLRRDSHKV DNFLKVLRCQ IAHQNNC (197)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	22.4 kDa
<b>Entrez</b>	<a href="#">19109</a>
<b>Expiration Date</b>	6 months from date of receipt.