

Product Datasheet

Canine CCL17 Recombinant Protein (orb1972940)

| | |
|-------------------------|---|
| Description | The Canine CCL17 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th... |
| Target | CCL17 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | ARGTNVGREC CLEYFKGAIP ISRLTRWYKT SGECPKDAIV FVTVQGKSIC SDPKDKRVKK AVRYLQRTWK GGPQES (76) |
| Purity | 98% |
| Source | Yeast |
| MW | 8.6 kDa |
| Entrez | 403586 |
| Expiration Date | 6 months from date of receipt. |