

Product Datasheet

Canine CCL17 Recombinant Protein (orb1972941)

Description	The Canine CCL17 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	CCL17
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ARGTNVGREC CLEYFKGAIP ISRLTRWYKT SGECPKDAIV FVTVQGKSIC SDPKDKRVKK AVRYLQRTWK GGPQES (76)
Purity	98%
Source	Yeast
MW	8.6 kDa
Entrez	403586
Expiration Date	6 months from date of receipt.