

## Product Datasheet

### Equine IL-36RA Recombinant Protein (orb1972960)

<b>Description</b>	The Equine IL-36 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturall...
<b>Target</b>	IL-36RA
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	VLSGALCFRM KDAALKVLYL HDNQLLAGGL HTGKVMKGEE ISVVPNRFLD AKLSPVILGV QGGSQLSCG TEQEPTLKLE PVNIMELYLG AKEYKSFTFY RRDTGLTSSF ESAAYPGWFL CTVPEADQPL RLTQLPEDAS GDDPITDFYF QQCD (154) [Asp137Glu]
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.0 kDa
<b>Entrez</b>	<a href="#">100065154</a>
<b>Expiration Date</b>	6 months from date of receipt.