

Product Datasheet

Equine IL-36RA Recombinant Protein (orb1972961)

Description	The Equine IL-36 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturall...
Target	IL-36RA
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VLSGALCFRM KDAALKVLYL HDNQLLAGGL HTGKVMKGEE ISVVPNRFLD AKLSPVILGV QGGSQLSCG TEQEPTLKLE PVNIMELYLG AKEYKSFTFY RRDGTGLTSSF ESAAYPGWFL CTVPEADQPL RLTQLPEDAS GDDPITDFYF QQCD (154) [Asp137Glu]
Purity	98%
Source	Yeast
MW	17.0 kDa
Entrez	100065154
Expiration Date	6 months from date of receipt.