

Product Datasheet

Ferret IL-6 Recombinant Protein (orb1972962)

Description	The Ferret IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-6
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	FPTPGPLGGD SKDDATSNRP PLTSADKMED FIRFILGKIS ALKKEMCEKY NKCEDSKAL AENNLNLPKL AEEDKCFQSQ FNQETCLTRI TTGLQEFQIH LKYLEANYEG NKNNAHSVYI STKHLLQKLR PMNRVEVTTTP DPTDSSLQA LFKSQDKWLK HVTIHLILRS LEDFLQFSLR AIRIM (185)
Purity	98%
Source	Yeast
MW	21.2 kDa
Entrez	101671465
Expiration Date	6 months from date of receipt.