

## Product Datasheet

# Cynomolgus Monkey TIM-3 Recombinant Protein (orb1972964)

<b>Description</b>	The Cynomolgus Monkey TIM-3 yeast-derived recombinant protein is not tagged and is naturally endotox...
<b>Target</b>	TIM-3
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	SEVEYIAEVG QNAYLPCSYT PAPPGNLVPV CWGKGACPVF DCSNVLRD NRDVNDRTSG RYWLKGDFFK GDVSLTIENV TLADSGVYCC RIQIPGIMND EKHNVKLVVI KPAKVTPAPT LQRDLTSAFP RMLTTGEHGP AETQTPGSLP DVNLTVSNNF CELQIFLTN ELRDSGAT (178)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	19.4 kDa
<b>Entrez</b>	Not Assigned
<b>Expiration Date</b>	6 months from date of receipt.