

Product Datasheet

Cynomolgus Monkey TIM-3 Recombinant Protein (orb1972966)

Description	The Cynomolgus Monkey TIM-3 yeast-derived recombinant protein is not tagged and is naturally endotox...
Target	TIM-3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SEVEYIAEVG QNAYLPCSYT PAPPGNLVPV CWGKGACPVF DCSNVVLRTD NRDVNDRTSG RYWLKGDFFK GDVSLTIENV TLADSGVYCC RIQIPGIMND EKHNVKLVVI KPAKVTPAPT LQRDLTSAFP RMLTTGEHGP AETQTPGSLP DVNLTVSNFF CELQIFTLTN ELRDSGAT (178)
Purity	98%
Source	Yeast
MW	19.4 kDa
Entrez	Not Assigned
Expiration Date	6 months from date of receipt.