

Product Datasheet

Dolphin TIM-3 Recombinant Protein (orb1973064)

Description	The Dolphin TIM-3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VYVSEVGQDA DLPCAYSPPT PEDLVPVCWG KGPCPLFQCH SLVLSTDGRN LSYWTSSRYL LKRNFRKGDV TLTIKNVTLE DSGTYCCRIQ FPGLMNDQKS NLELVIKPAK VTPAGTPWRD ITAAFPRMLT TEGPGSETQT LETLHDKNQT EISTLANELQ DMGATTRI (168)
Purity	98%
Source	Yeast
MW	18.6 kDa
Entrez	101331594
Expiration Date	6 months from date of receipt.