

Product Datasheet

Rat Prolactin Receptor Recombinant Protein (orb1973066)

Description	The Rat Prolactin Receptor yeast-derived recombinant protein is not tagged and is naturally endotoxi...
Target	PRLR
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QSPPGKPEIH KCRSPDKETF TCWWNPGTDG GLPTNYSLTY SKEGEKTTYE CPDYKTS GPN SCFFSKQYTS IWKIYIITVN ATNQMGS SSS DPLYVDVTYI VEPEPPRNL T LEVKQLKDKK TYLWVKWSPP TITDVKTGWF TMEYEIRLKP EEAEWEIHF TGHQTQFKVF DLYPGQKYL V QTRCKPDHGY WSRWSQESSV EMPNDFTLKD (210)
Purity	98%
Source	Yeast
MW	24.6 kDa
Entrez	24684
Expiration Date	6 months from date of receipt.