

Product Datasheet

Rat Prolactin Receptor Recombinant Protein (orb1973067)

Description	The Rat Prolactin Receptor yeast-derived recombinant protein is not tagged and is naturally endotoxi...
Target	PRLR
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QSPPGKPEIH KCRSPDKETF TCWWNPGTDG GLPTNYSLTY SKEGEKTTYE CPDYKTSGPN SCFFSKQYTS IWKIYITVN ATNQMGSSSS DPLYVDVTYI VEPEPPRNLTV LEVKQLKDKK TYLWVKWSPP TITDVKTGWF TMEYEIRLKP EEAEWEIHF TGHQTQFKVF DLYPGQKYL QTRCKPDHGY WSRWSQESSV EMPNDFTLKD (210)
Purity	98%
Source	Yeast
MW	24.6 kDa
Entrez	24684
Expiration Date	6 months from date of receipt.