

## Product Datasheet

# Cynomolgus Monkey Fas Ligand Recombinant Protein (orb1973069)

<b>Description</b>	The Cynomolgus Monkey Fas Ligand yeast-derived recombinant protein is not tagged and is naturally en...
<b>Target</b>	Fas Ligand
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	QLFHLQKELA ELRESTSQKH TASSLEKQIG HPSPPPEKKE QRKVAHLTGK PNSRSMPLEW EDTYGIVLLS GVKYKKGGLV INETGLYFVY SKVYFRGQSC TNLPLSHKVY MRNSKYPQDL VMMEGKMMSY CTTGQMWAHS SYLGAVFNLT SADHLYVNVS ELSLVNFEES QTFGLYKL (179)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	20.4 kDa
<b>Entrez</b>	<a href="#">102139406</a>
<b>Expiration Date</b>	6 months from date of receipt.