

Product Datasheet

Swine/Bat Pleiotrophin Recombinant Protein (orb1973075)

Description	The Swine/Bat Pleiotrophin yeast-derived recombinant protein is not tagged and is naturally endotoxi...
Target	Pleiotrophin
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GKKEKPEKKV KKSDCGEWQW SVCVPTSGDC GLGTREGTRT GAECKQTMKT QRCKIPC�WK KQFGAECKYQ FQAWGECDLN TALKTRTGSL KRALHNADCQ KTVTISKPCG KVTKPKPQAE SKKKKKEGKK QEKMLD (136)
Purity	98%
Source	Yeast
MW	15.3 kDa
Entrez	397609
Expiration Date	6 months from date of receipt.