

Product Datasheet

Feline IL-13 Recombinant Protein (orb1973080)

| | |
|-------------------------|--|
| Description | The Feline IL-13 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th... |
| Target | IL-13 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | SPGPHSRREL KELIEELVNI TQNQVSLCNG SMVWSVNLTT GMQYCAALES LINVSDCTAI QRTQRMLKAL CTQKPSAGQT ASERSRDTKI EVIQLVKNLL NHLRRNFRHG NFK (113) |
| Purity | 98% |
| Source | Yeast |
| MW | 12.7 kDa |
| Entrez | 101084678 |
| Expiration Date | 6 months from date of receipt. |