

Product Datasheet

Feline IL-13 Recombinant Protein (orb1973082)

Description	The Feline IL-13 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	IL-13
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPGPHSRREL KELIEELVNI TQNQVSLCNG SMVWSVNLTT GMQYCAALES LINVSDCTAI QRTQRMLKAL CTQKPSAGQT ASERSRDTKI EVIQLVKNLL NHLRRNFRHG NFK (113)
Purity	98%
Source	Yeast
MW	12.7 kDa
Entrez	101084678
Expiration Date	6 months from date of receipt.