

Product Datasheet

Mink IL-8 Recombinant Protein (orb1973083)

Description	The Mink IL-8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The p...
Target	IL-8
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AVLSRVSSEL RCQCIKTHST PFHPKYIKEL RVIDSGPHCE NSEIIVKLVN GKEVCLDPKE DWVQRIVQIF LKKA EKQAS (79)
Purity	98%
Source	Yeast
MW	9.0 kDa
Entrez	448792
Expiration Date	6 months from date of receipt.