

Product Datasheet

Mink IL-8 Recombinant Protein (orb1973085)

| | |
|-------------------------|---|
| Description | The Mink IL-8 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The p... |
| Target | IL-8 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | AVLSRVSSEL RCQCIKTHST PFHPKYIKEL RVIDSGPHCE NSEIIVKLVN GKEVCLDPKE DWVQRIVQIF LKKAQEKQAS (79) |
| Purity | 98% |
| Source | Yeast |
| MW | 9.0 kDa |
| Entrez | 448792 |
| Expiration Date | 6 months from date of receipt. |