



Product Datasheet

Factor IX / F9 antibody [PB-3] (orb1973109)


Description	Factor IX / F9 antibody [PB-3]
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Target	Factor IX / F9 antibody [PB-3]
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C Phosphate buffered solution, pH 7.2, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Clonality	Monoclonal
Dilution Range	Western Blot 1:500-1:1000 Immunofluorescence 1:200 - 1:300 Immunohistochemistry (Paraffin) 1:200 - 1:300
Expiration Date	12 months from date of receipt.



This blot was produced using a 12% SDS-P...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...