



Product Datasheet

Influenza A Virus Nucleoprotein, Detection antibody [BS772] (orb1973129)

Description	Influenza A Virus Nucleoprotein, Detection antibody [BS772]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, LF, WB
Target	Influenza A Virus Nucleoprotein, Detection antibody [BS772]
Storage	Phosphate-buffered solutions, pH 7.2. Long time storage is recommended at -20°C Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.
Note	For research use only
Clonality	Monoclonal
Dilution Range	Western Blot 1:2000-1:4000 ELISA 1:5000 - 1:10000 Lateral flow
Expiration Date	12 months from date of receipt.



Detection of Influenza A Virus and Influenza...



Detection of Influenza A Virus and Influenza...