



## Product Datasheet

# Influenza A Virus Nucleoprotein, Capture antibody [BS771] (orb1973130)

<b>Description</b>	Influenza A Virus Nucleoprotein, Capture antibody [BS771]
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, LF, WB
<b>Target</b>	Influenza A Virus Nucleoprotein, Capture antibody [BS771]
<b>Storage</b>	Phosphate-buffered solutions, pH 7.2. Long time storage is recommended at -20°C Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	Western Blot 1:2000-1:5000 ELISA 1:5000 - 1:10000 Lateral flow
<b>Expiration Date</b>	12 months from date of receipt.



Detection of Influenza A Virus and Influenza...



Detection of Influenza A Virus and Influenza...