


Product Datasheet

SARS-CoV-2 (COVID-19) Spike protein S2 Antibody (orb1973131)

Description	SARS-CoV-2 (COVID-19) Spike protein S2 Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IF, NeA, WB
Target	SARS-CoV-2 (COVID-19) Spike protein S2
Storage	Long time storage is recommended at -20°C PBS (Ph7.4), 20% Glycerol. 0.025% ProClin 300 was added as a preservative Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user
Note	For research use only
Clonality	Polyclonal
Dilution Range	Western Blot 1:2000 – 1:14000 Immunofluorescence 1:500 – 1:1000 Neutralizing Reported in the literature, not verified in house.
Expiration Date	12 months from date of receipt.



Western blot analysis of SARS-CoV-2 (COV...