


Product Datasheet

SARS-CoV-2 (COVID-19) Spike RBD antibody [10F] (orb1973133)

Description	SARS-CoV-2 (COVID-19) Spike RBD antibody [10F]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IF, NeA, WB
Target	SARS-CoV-2 (COVID-19) Spike RBD antibody [10F]
Storage	Long time storage is recommended at -20°C PBS (Ph7.4), 20% Glycerol. 0.025% ProClin 300 was added as a preservative Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user
Note	For research use only
Clonality	Monoclonal
Dilution Range	Western Blot 1:1000 - 1:12000 Immunofluorescence 1:500 - 1:1000 Neutralizing Reported in the literature, not verified in house.
Expiration Date	12 months from date of receipt.



Western blot analysis of SARS-CoV-2 (COV...