


Product Datasheet

SARS-CoV-2 nucleocapsid protein, monoclonal [6C2D] Antibody (orb1973134)

Description	SARS-CoV-2 nucleocapsid protein, monoclonal [6C2D] Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Target	SARS-CoV-2 nucleocapsid protein, monoclonal [6C2D]
Storage	Long time storage is recommended at -20°C Phosphate-buffered solutions, pH 7.2. 50% Glycerol, 0.025% ProClin 300 was added as a preservative Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.
Note	For research use only
Clonality	Monoclonal
Dilution Range	Western Blot 1:2000-1:4000 ELISA 1:5000 - 1:10000
Expiration Date	12 months from date of receipt.



Nitrocellulose membrane was then blocked...