


Product Datasheet

SARS-CoV-2 nucleocapsid protein, Capture antibody [BS41] (orb1973136)

Description	SARS-CoV-2 nucleocapsid protein, Capture antibody [BS41]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, LF, WB
Target	SARS-CoV-2 nucleocapsid protein, Capture antibody [BS41]
Storage	Phosphate-buffered solutions, pH 7.2. Long time storage is recommended at -20°C Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.
Note	For research use only
Clonality	Monoclonal
Dilution Range	Western Blot 1:2000-1:4000 ELISA 1:5000 - 1:10000 Lateral flow
Expiration Date	12 months from date of receipt.



Detection of SARS-CoV-2 nucleocapsid pro...