



Product Datasheet

F4/80 [BM8], monoclonal Antibody (orb1973146)

Description	F4/80 [BM8], monoclonal Antibody
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Target	F4/80 [BM8], monoclonal
Storage	Store at -20°C Phosphate-buffered solution, pH 7.2, containing 0.025% ProClin 300, 20% Glycerol Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.
Note	For research use only
Clonality	Monoclonal
Dilution Range	Immunohistochemistry 1:100 Flow Cytometry 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Mouse peritoneal macrophages stained with...