



Product Datasheet

Carbonic Anhydrase IX / CA-IX antibody (orb1973156)


Description	Carbonic Anhydrase IX / CA-IX antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Synthetic peptide
Target	Carbonic Anhydrase IX / CA-IX antibody
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C 100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative. Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.
Note	For research use only
Clonality	Polyclonal
Dilution Range	Immunofluorescence 1:100 - 1:300 Immunohistochemistry (Paraffin) 1:50 - 1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...