


Product Datasheet

HA tag antibody [PB-4] (orb1973165)

Description	HA tag antibody [PB-4]
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Target	HA tag antibody [PB-4]
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C Phosphate buffered solution, pH 7.2, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Note	For research use only
Clonality	Monoclonal
Dilution Range	Western Blot 1:3000-1:5000 Immunofluorescence 1:500 - 1:1000
Expiration Date	12 months from date of receipt.



This blot was produced using a 7.5% SDS...