



Product Datasheet

GFP , monoclonal antibody (orb1973166)

Description	GFP , monoclonal antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Target	GFP , monoclonal antibody
Storage	Long time storage is recommended at -20°C Phosphate buffered solution, pH 7.2, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.
Note	For research use only
Clonality	Monoclonal
Dilution Range	Western Blot 1:1000 - 1:2000 Immunofluorescence 1:300 - 1:500
Expiration Date	12 months from date of receipt.



GFP, monoclonal antibody at 1/1000 dilut...



Immunofluorescence staining of methanol...