


## Product Datasheet

### GST tag antibody [PA-3] (orb1973167)

<b>Description</b>	GST tag antibody [PA-3]
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IP, WB
<b>Target</b>	GST tag antibody [PA-3]
<b>Storage</b>	Store at +4°C for short term storage. Long time storage is recommended at -20°C Phosphate buffered solution, pH 7.2, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	Western Blot 1:2000-1:5000 Immunofluorescence 1:500 - 1:1000
<b>Expiration Date</b>	12 months from date of receipt.



This blot was produced using a 10% SDS-P...