

Product Datasheet

CD38 inhibitor 78c (orb1973207)

Description	CD38 inhibitor 78c
Form/Appearance	Dry powder
Storage	Store in dry, dark place for 12 months.
Note	For research use only
Formula	$C_{22}H_{27}N_3O_3S$
Purity	≥98%
MW	413.53 kDa
CAS Number	1700637-55-3
SMILES	<chem>COCCO[C@@H]1CC[C@H](CC1)NC1=CC(=O)N(C)C2C=CC(=CC=21)C1=CN=CS1</chem>
Solubility (25°C)	Soluble in DMSO
Expiration Date	12 months from date of receipt.