

Product Datasheet

Filastatin (orb1973340)

| | |
|--------------------------|---|
| Description | Filastatin |
| Form/Appearance | Dry powder |
| Storage | Store in dry, dark place for 12 months. |
| Note | For research use only |
| Formula | $C_{18}H_{18}ClN_3O_3$ |
| Purity | ≥98% |
| MW | 359.81 kDa |
| CAS Number | 431996-53-1 |
| SMILES | <chem>CC1C=CC(=CC=1Cl)C(=O)N1CCN(CC1)C1C=CC(=CC=1)[N+][O-]=O</chem> |
| Solubility (25°C) | Soluble in DMSO |
| Expiration Date | 12 months from date of receipt. |