

Product Datasheet

Benclonidine (orb1973421)

Description	Benclonidine
Form/Appearance	Dry powder
Storage	Store in dry, dark place for 12 months.
Note	For research use only
Formula	$C_{16}H_{13}Cl_2N_3O$
Purity	$\geq 98\%$
MW	334.2 kDa
CAS Number	57647-79-7
SMILES	<chem>O=C(C1C=CC=CC=1)N1CCN=C1NC1C(Cl)=CC=CC=1Cl</chem>
Solubility (25°C)	Soluble in DMSO
Expiration Date	12 months from date of receipt.