

Product Datasheet

Pirfenidone D5 (orb1973508)

Description	Pirfenidone D5
Form/Appearance	Dry powder
Storage	Store in dry, dark place for 12 months.
Note	For research use only
Formula	$C_{12}H_{11}NO$
Purity	$\geq 98\%$
MW	190.25 kDa
CAS Number	1020719-62-3
SMILES	<chem>[2H]C1=C(C([2H])=C([2H])C([2H])=C1[2H])N1C=C(C)C=CC1=O</chem>
Solubility (25°C)	Soluble in DMSO
Expiration Date	12 months from date of receipt.