

Product Datasheet

Pirfenidone (orb1973521)

Description	Pirfenidone
Form/Appearance	Dry powder
Storage	Store in dry, dark place for 12 months.
Note	For research use only
Formula	$C_{12}H_{11}NO$
Purity	$\geq 98\%$
MW	185.22 kDa
CAS Number	53179-13-8
SMILES	<chem>CC1=CN(C(=O)C=C1)C1C=CC=CC=1</chem>
Solubility (25°C)	Soluble in DMSO, not in water
Expiration Date	12 months from date of receipt.