

Product Datasheet

TAK-243 (MLN7243) (orb1973963)

Description	TAK-243 (MLN7243)
Form/Appearance	Dry powder
Storage	Store in dry, dark place for 12 months.
Note	For research use only
Formula	$C_{19}H_{20}F_3N_5O_5S_2$
Purity	≥98%
MW	519.52 kDa
CAS Number	1450833-55-2
SMILES	<chem>NS(=O)(=O)OC[C@H]1C[C@@H](NC2=CC=NC3=CC(=NN3)C2=CC(=C=C2)SC(F)(F)F)[C@H](O)[C@@H]1O</chem>
Solubility (25°C)	DMSO up to 100 mM
Expiration Date	12 months from date of receipt.