

Product Datasheet

Recombinant human C1Q protein, N-His (HEK293) (orb1974300)

Description	Recombinant human C1Q protein, N-His (HEK293)
Endotoxins	Not analyzed
Target	C1QA
Form/Appearance	Liquid
Concentration	> 1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	N-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	> 90% as determined by SDS-PAGE
MW	48.6 kDa
Uniprot ID	P02747
Expiration Date	6 months from date of receipt.