

Product Datasheet

Recombinant human Coagulation factor II protein/F2, C-His (HEK293) (orb1974301)

Description	Recombinant human Coagulation factor II protein/F2, C-His (HEK293)
Endotoxins	Not analyzed
Target	F2
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	44-413/622
Purity	> 90% as determined by SDS-PAGE
MW	31.8 kDa
Uniprot ID	P00734
Expiration Date	6 months from date of receipt.