

Product Datasheet

Recombinant human FCER1A protein, C-His-Avi (HEK293) (orb1974306)

| | |
|--------------------------|--|
| Description | Recombinant human FCER1A protein, C-His-Avi (HEK293) |
| Endotoxins | Not analyzed |
| Target | FCER1A |
| Form/Appearance | Liquid |
| Concentration | > 1mg/ml |
| Storage | Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles. |
| Tag | C-His-Avi |
| Note | For research use only |
| Application notes | Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water |
| Purity | > 90% as determined by SDS-PAGE |
| MW | 1.6 kDa |
| Uniprot ID | P12319 |
| Expiration Date | 6 months from date of receipt. |