

Product Datasheet

Recombinant human FCER2/CD23 protein,N-His-Avi (HEK293) (orb1974307)

Description	Recombinant human FCER2/CD23 protein,N-His-Avi (HEK293)
Endotoxins	Not analyzed
Target	FCER2
Form/Appearance	Liquid
Concentration	> 1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	N-His-Avi
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	> 90% as determined by SDS-PAGE
MW	34.6 kDa
Uniprot ID	P06734
Expiration Date	6 months from date of receipt.