

Product Datasheet

Recombinant human Angiopoietin-2 protein, N-His-Avi (HEK293) (orb1974314)

Description	Recombinant human Angiopoietin-2 protein, N-His-Avi (HEK293)
Endotoxins	< 1.0 EU/μg as determined by LAL
Target	ANGPT2
Form/Appearance	Lyophilized or liquid
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	N-His-Avi tag
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	> 95% as determined by Tris-Bis PAGE; > 95% as determined by SEC-HPLC
Uniprot ID	O15123
Expiration Date	6 months from date of receipt.