

Product Datasheet

Recombinant mouse FCGR1A/CD64 protein, (HEK293) (orb1974338)

Description	Recombinant mouse FCGR1A/CD64 protein, (HEK293)
Endotoxins	Not analyzed
Target	FCGR1A
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His-Avi
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	25-297/404
Purity	> 90% as determined by SDS-PAGE
MW	34 kDa
Uniprot ID	P26151
Expiration Date	6 months from date of receipt.