

Product Datasheet

Recombinant human Fc gamma RIIIB / CD16b (NA2) Protein; C-His (HEK293) (orb1974339)

Description	Recombinant human Fc gamma RIIIB / CD16b (NA2) Protein; C-His (HEK293)
Endotoxins	Not analyzed
Target	FCGR3B
Form/Appearance	Lyophilized or liquid
Concentration	> 1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	> 90% as determined by SDS-PAGE
MW	23.2 kDa
Uniprot ID	O75015
Expiration Date	6 months from date of receipt.