

Product Datasheet

Recombinant human Fc gamma RIIIA/CD16a (F176) Protein, C-His(HEK293) (orb1974346)

Description	Recombinant human Fc gamma RIIIA/CD16a (F176) Protein, C-His(HEK293)
Endotoxins	Not analyzed
Target	FCGR3A
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	> 90% as determined by SDS-PAGE
MW	24 kDa
Uniprot ID	P08637
Expiration Date	6 months from date of receipt.