

Product Datasheet

Recombinant human Vasopressin-neurophysin 2-copeptin protein (HEK293) (orb1974348)

Description	Recombinant human Vasopressin-neurophysin 2-copeptin protein (HEK293)
Endotoxins	Not analyzed
Target	AVP
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	20-164/164
Purity	> 90% as determined by SDS-PAGE
MW	17.4 kDa
Uniprot ID	P01185
Expiration Date	6 months from date of receipt.