

Product Datasheet

Recombinant human SHBG protein (HEK293), His (orb1974360)

Description	Recombinant human SHBG protein (HEK293), His
Endotoxins	Not analyzed
Target	SHBG
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	30-402/402
Purity	> 90% as determined by SDS-PAGE
MW	41.5 kDa
Uniprot ID	A0A6P7EHR7
Expiration Date	6 months from date of receipt.