

Product Datasheet

Recombinant Human Fc gamma RIIB/C(CD32b/c) (HEK293), His-Avi (orb1974361)

Description	Recombinant Human Fc gamma RIIB/C(CD32b/c) (HEK293), His-Avi
Endotoxins	Not analyzed
Target	FCGR2B
Form/Appearance	Liquid
Concentration	> 0.5mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	> 90% as determined by SDS-PAGE
MW	21.4 kDa
Uniprot ID	P31994
Expiration Date	6 months from date of receipt.